

Anti-HPV-16 E6 antibody [HPV-13E2] (STJ16101056)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

