

## Anti-HPV-18 E6 antibody [HPV-4G3] (STJ16101058)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IHC          |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.  |
| Conjugation         | Unconjugated  |
| Purification        | Affinity purified from tissue culture.  |
| Dilution Range      | ELISA (solid phase: not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date). |
| Formulation         | PBS with 0.02% Sodium Azide.  |
| Isotype             | IgG1k   |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt.   |

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

