

Anti-MFG-E8/Lactadherin/p47 antibody [MFG-06] (STJ16101070)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/IF/WB
Host / Source	Mouse
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide
Isotype	IgG1K
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4240
Gene Symbol	MFG8
UniProt ID	MFGM_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

