

Anti-Cdk1/p34cdc2 antibody [POH-1] (STJ16101079)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB/IF/IHC
Host / Source	Mouse
Reactivity	Cow/Human/Mink/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow Cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunoblotting (1-2 µg/ml). Immunofluorescence (1-2 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min at RT; requires boiling tissue sections in 10mM Citrate Buffer, pH 6.0, for 10-2
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2ak
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	983
Gene Symbol	CDK1
UniProt ID	CDK1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

