

Anti-MAP3K1 antibody [EBS-T-008] (STJ16101082)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/WB/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow Cytometry (0, 5-1 µg/million cells in 0.1 ml). Immunoblotting (0, 5-1 µg/ml). Immunohistology (0, 5-1, 0 µg/ml for 30 minutes at RT; Staining of formalin-fixed tissue
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2ak
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	4214
Gene Symbol	MAP3K1
UniProt ID	M3K1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

