

## Anti-CD25/IL2RA antibody [EBS-CD-017] (STJ16101100)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: 0, 1-100 µg/ml; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow Cytometry (1-2 Mu g/million cells in 0, 1 ml for 30 min, at 4°C). Immunofluorescence (1-2 Mu g/ml). Immunohistology (1-2 Mu g/ml for 30 min at RT; no suitable antigen)
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2ak
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---