

## Anti-CD63 antibody [EBS-CD-036] (STJ16101120)

| GENERAL INFORMATION    |                                                                                                                                                                                                                                       |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type           | Primary antibodies                                                                                                                                                                                                                    |
| Applications           | FC/IF/IHC                                                                                                                                                                                                                             |
| Host / Source          | Mouse                                                                                                                                                                                                                                 |
| Reactivity             | Human/Mouse                                                                                                                                                                                                                           |
| PRODUCT PROPERTIES     |                                                                                                                                                                                                                                       |
| Clonality              | Monoclonal                                                                                                                                                                                                                            |
| Concentration          | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.                                                                                                                                                                                    |
| Conjugation            | Unconjugated                                                                                                                                                                                                                          |
| Purification           | Affinity purified from tissue culture.                                                                                                                                                                                                |
| Dilution Range         | Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (0, 5-1, 0 µg/ml for 30 min at RT; no special pre-treatment is required for IHC staining of formalin/paraffin tissues). |
| Formulation            | PBS with 0.02% Sodium Azide.                                                                                                                                                                                                          |
| Isotype                | IgG1k                                                                                                                                                                                                                                 |
| Storage Instruction    | Store for up to 1 year at 2-8°C upon receipt.                                                                                                                                                                                         |
| TARGET INFORMATION     |                                                                                                                                                                                                                                       |
| Gene ID                | <a href="#">967</a>                                                                                                                                                                                                                   |
| Gene Symbol            | <a href="#">CD63</a>                                                                                                                                                                                                                  |
| UniProt ID             | <a href="#">CD63_HUMAN</a>                                                                                                                                                                                                            |
| ADDITIONAL INFORMATION |                                                                                                                                                                                                                                       |
| Note                   | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.                                                                                                               |