

Anti-CD66e-CEACAM5 antibody [CB-30] (STJ16101121)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/FC/IF/IHC
Host / Source	Mouse
Reactivity	Human/Monkey
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	ELISA (solid phase: not known; tracer: 0, 001-100 µg/ml for 30 min at RT). Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues r
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG1k
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.
TARGET INFORMATION	
Gene ID	1048
Gene Symbol	CEACAM5
UniProt ID	CEAM5_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.