

Anti-CD80 antibody [C80-99] (STJ16101125)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IHC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30-60 min at RT; for staining of formalin-fixed tissues no suitable antigen retrieval method is known to date).
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2bk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.

TARGET INFORMATION

Gene ID	941
Gene Symbol	CD80
UniProt ID	CD80_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---