

Anti-Interleukin-10 antibody (STJ16101143)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/ELISA/IF/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Affinity purified
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	3586
Gene Symbol	IL10
UniProt ID	IL10_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---