

Anti-Desmin K5 antibody (STJ16101153)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB
Host / Source	Rabbit
Reactivity	Human/Bovine/Caprine/Chicken/Hamster/Mouse/Rat/Zebrafish

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	IHC-F 1:25-1:100 IHC-P 1:25-1:100 WB 1:100-1:500
Formulation	PBS with 0.09% sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

