

Anti-TNF-Alpha antibody (STJ16101171)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/ELISA/IHC-P/WB
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	100 Mu g/mL
Conjugation	Unconjugated
Purification	Affinity purified
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	Ig
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	7124
Gene Symbol	TNF
UniProt ID	TNFA_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---