

Anti-caspase-1 antibody (C-Term) (STJ16101311)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P
Host / Source	Rabbit

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	1.0 mg/mL
Conjugation	Unconjugated
Purification	Affinity purified
Formulation	PBS with 0.02% sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	834
Gene Symbol	CASP1
UniProt ID	CASP1_HUMAN
Immunogen Region	C-Term

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---