

## Anti-HSV1 antibody (STJ16101333)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                   |
|---------------------|-------------------------------------------------------------------|
| Clonality           | Polyclonal                                                        |
| Conjugation         | Unconjugated                                                      |
| Purification        | Affinity purified                                                 |
| Dilution Range      | IHC 1:100-1:500                                                   |
| Formulation         | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt.                    |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

