

## Anti-E-Cadherin antibody [BS38] (STJ16101355)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human/Dog/Rat/Mouse/Pig/Sheep/Rabbit

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:100-1:400
Formulation	TRIS with 0.03% sodium azide, pH7, 2
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">999</a>
Gene Symbol	<a href="#">CDH1</a>
UniProt ID	<a href="#">CADH1_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

