

Anti-Pan Cytokeratin antibody [BS5] (STJ16101360)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F
Host / Source	Mouse
Reactivity	Human/Sheep/Pig/Dog/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:100-1:400 IHC-F
Formulation	TRIS with 0.03% sodium azide, pH7, 2
Isotype	IgG1
Storage Instruction	Store at 2-8°C upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---