

Anti-UP II antibody (STJ16101373)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Rabbit, Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Mix
Concentration	Lot specific
Conjugation	Unconjugated
Dilution Range	IHC-P 1:50-1:100
Formulation	Tissue culture supernatant with sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

