

Anti-Cytokeratin 7, Keratin K7 antibody [RCK105] (STJ16101384)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/IHC-F/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Zebrafish/Dog/Caprine/Feline/Hamster

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Formulation	PBS with 0.09% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

