

## Mouse IgG2a Isotype Control [MOPC-173] {Biotin} (STJ16101393)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/ICC/WB/IP/IHC-P/IHC-F/ELISA
Host / Source	Mouse

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Biotin
Formulation	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---