

Anti-ABCB1 antibody [JSB-1] (STJ16101408)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Chinese Hamster

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	100 µg/mL
Conjugation	Unconjugated
Dilution Range	FC 1:10 IHC-F IHC-P 1:20 WB
Formulation	PBS with 0.1% BSA and 0.02% sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	5243
Gene Symbol	ABCB1
UniProt ID	MDR1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

