

Anti-PAX-2 antibody [EP235] (STJ16101670)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P |
| Host / Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Dilution Range | IHC-P 1:10-1:50 |
| Formulation | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

