

Anti-CA19-9 antibody [C241:5:1:4] (STJ16101686)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:200 IHC-F
Formulation	Purified in PBS, 1% BSA and 15 mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

