

Anti-Calcitonin antibody [CL1948] (STJ16101688)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	>=29 mg/L
Conjugation	Unconjugated
Dilution Range	IHC-P 1:200
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG2b
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---

