

Anti-Caveolin-1 antibody [4D6] (STJ16101692)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:200 WB
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG2ak
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	857
Gene Symbol	CAV1
UniProt ID	CAV1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

