

Anti-Delta-Sarcoglycan antibody [omegaSarc3/12C1] (STJ16101728)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit/Dog/Chicken/Hamster/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Lyophilized
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	6444
Gene Symbol	SGCD
UniProt ID	SGCD_HUMAN

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

