

Anti-Dystrophin antibody (Rod Domain) [DY4/6D3] (STJ16101729)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Mouse
Reactivity	Human/Mouse/Rat/Rabbit/Dog/Hamster/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Lyophilized
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	1756
Gene Symbol	DMD
UniProt ID	DMD_HUMAN
Immunogen Region	Rod Domain

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

