

## Anti-Dystrophin antibody (C-Term) [DY8/6C5] (STJ16101730)

### GENERAL INFORMATION

|               |  |
|---------------|--|
| Product Type  | Primary antibodies                         |
| Applications  | IHC-F                                      |
| Host / Source | Mouse                                      |
| Reactivity    | Human/Mouse/Rat/Rabbit/Dog/Hamster/Chicken |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal                                     |
| Conjugation         | Unconjugated                                   |
| Formulation         | Lyophilized                                    |
| Isotype             | IgG1   |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### TARGET INFORMATION

|                  |                           |
|------------------|---------------------------|
| Gene ID          | <a href="#">1756</a>      |
| Gene Symbol      | <a href="#">DMD</a>       |
| UniProt ID       | <a href="#">DMD_HUMAN</a> |
| Immunogen Region | C-Term                    |

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

