

Anti-Dystrophin antibody (N-Term) [DY10/12B2] (STJ16101731)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-F |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Formulation | Lyophilized |
| Isotype | IgG2a |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

TARGET INFORMATION

| | |
|------------------|---------------------------|
| Gene ID | 1756 |
| Gene Symbol | DMD |
| UniProt ID | DMD_HUMAN |
| Immunogen Region | N-Term |

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

