

## Anti-CD34 antibody [QBEnd/10] (STJ16101737)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | >=261 mg/L                                     |
| Conjugation         | Unconjugated                                   |
| Dilution Range      | IHC-P 1:100                                    |
| Formulation         | Tissue culture supernatant with Sodium azide   |
| Isotype             | IgG1                                           |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">947</a>        |
| Gene Symbol | <a href="#">CD34</a>       |
| UniProt ID  | <a href="#">CD34_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

