

Anti-Fascin antibody [IM20] (STJ16101740)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:200-1:400 WB 1:100-1:200
Formulation	Tissue culture supernatant with 15 mM Sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	6624
Gene Symbol	FSCN1
UniProt ID	FSCN1_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

