

## Anti-Granzyme B antibody [11F1] (STJ16101744)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	>=47 mg/L
Conjugation	Unconjugated
Dilution Range	IHC-P 1:80
Formulation	Tissue culture supernatant with Sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">3002</a>
Gene Symbol	<a href="#">GZMB</a>
UniProt ID	<a href="#">GRAB_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

