

Anti-Merosin Laminin Alpha 2 Chain antibody [Mer3/22B2] (STJ16101762)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F
Host / Source	Mouse
Reactivity	Human/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Formulation	Lyophilized
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---