

## Anti-MITF antibody [34CA5] (STJ16101768)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:10-1:20 IHC-F
Formulation	Tissue culture supernatant with 15mM Sodium azide
Isotype	IgG1k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">4286</a>
Gene Symbol	<a href="#">MITF</a>
UniProt ID	<a href="#">MITF_HUMAN</a>

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

