

## Anti-MSH6 antibody [PU29] (STJ16101771)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	>=200 mg/L
Conjugation	Unconjugated
Dilution Range	IHC-P 1:100
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">2956</a>
Gene Symbol	<a href="#">MSH6</a>
UniProt ID	<a href="#">MSH6_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

