

Anti-Prostatic Acid Phosphatase antibody [PASE4LJ] (STJ16101781)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:100 IHC-F
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	55
Gene Symbol	ACP3
UniProt ID	PPAP_HUMAN

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

