

## Anti-Retinoblastoma Gene Protein antibody [13A10] (STJ16101790)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/IHC-F/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Dilution Range	IHC-P 1:25-1:50 IHC-F WB 1:250
Formulation	Tissue culture supernatant with 15mM sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

