

Anti-Thyroid Transcription Factor-1 antibody [SPT24] (STJ16101797)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IHC-P |
| Host / Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|--|
| Clonality | Monoclonal |
| Concentration | >=108 mg/L |
| Conjugation | Unconjugated |
| Dilution Range | IHC-P 1:200 |
| Formulation | Tissue culture supernatant with sodium azide |
| Isotype | IgG1 |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.

