

## Anti-Vimentin antibody [SRL33] (STJ16101801)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	>=60 mg/L
Conjugation	Unconjugated
Dilution Range	IHC-P 1:400
Formulation	Tissue culture supernatant with 15 mM sodium azide
Isotype	IgG2a
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### TARGET INFORMATION

Gene ID	<a href="#">7431</a>
Gene Symbol	<a href="#">VIM</a>
UniProt ID	<a href="#">VIME_HUMAN</a>

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**

