

## Anti-von Willebrand Factor antibody [36B11] (STJ16101802)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IHC-P              |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                |
|---------------------|------------------------------------------------|
| Clonality           | Monoclonal                                     |
| Concentration       | >=21 mg/L                                      |
| Conjugation         | Unconjugated                                   |
| Dilution Range      | IHC-P 1:200                                    |
| Formulation         | Tissue culture supernatant with sodium azide   |
| Isotype             | IgG2a                                          |
| Storage Instruction | Store at 2-8°C for up to 1-year, upon receipt. |

### TARGET INFORMATION

|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">7450</a>      |
| Gene Symbol | <a href="#">VWF</a>       |
| UniProt ID  | <a href="#">VWF_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|

