

## Anti-CM5 antibody (STJ16101825)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/WB
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity purified
Dilution Range	IHC-P 1:200 WB 1:2000-1:4000
Formulation	PBS with 1% BSA and sodium azide
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---