

Anti-CD221 IGF Receptor antibody [1H7] (STJ16101843)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/FC/WB
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified
Dilution Range	IHC-P 1:10-1:40 FC 1:10 WB
Formulation	PBS with azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---