

Anti-MART-1 (MELAN A) antibody [M2-9E3] (STJ16101844)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC-F/IHC-P
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-F IHC-P 1:100-1:200
Formulation	Bioreactor Concentrate with 0.05% Azide
Isotype	IgG2bk
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---