

Anti-LMO2 antibody [1A9-1] (STJ16101845)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	200 Mu g IgG1/mL
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:100-1:200
Formulation	PBS with < 0.1% sodium azide and 0.1% gelatin
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4005
Gene Symbol	LMO2
UniProt ID	RBTN2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---