

Anti-MDM2 antibody [SMP14] (STJ16101854)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/IHC-P/WB
Host / Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified
Dilution Range	IHC-F IHC-P 1:50-1:100 WB
Formulation	Affinity purified antibody in aqueous buffered solution containing $\leq 0.09\%$ sodium azide
Isotype	IgG1
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	4193
Gene Symbol	MDM2
UniProt ID	MDM2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---