

## Anti-CD43 antibody [DF-1] (STJ16101856)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:100-1:200 WB 1:100-1:200
Formulation	Bioreactor Concentrate with 0.05% Azide
Isotype	IgG2k
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---