

Anti-E-CADHERIN antibody [MX020] (STJ16101857)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Concentrate
Dilution Range	IHC-P 1:100-1:200
Formulation	Tissue culture supernatant containing 15mM sodium azide
Isotype	IgG
Storage Instruction	Store at 2-8°C for up to 1-year, upon receipt.

TARGET INFORMATION

Gene ID	999
Gene Symbol	CDH1
UniProt ID	CADH1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---