

## Anti-CD45R0 antibody [EBS-CD-027] (STJ16101874)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IF/IHC          |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Concentration       | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.   |
| Conjugation         | Unconjugated   |
| Purification        | Affinity purified from tissue culture.   |
| Dilution Range      | Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (1-2 µg/ml for 30 min at RT; an appropriate antigen retrieval method for staining of formalin-fixed tissues has not been established to date |
| Formulation         | PBS with 0.02% Sodium Azide.   |
| Isotype             | IgG2ak   |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt.  |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|