

## Anti-CD45RA antibody [158-4D3] (STJ16101876)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IF/IHC          |
| Host / Source | Mouse              |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                                                                                                           |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                                                                                                                                                                                |
| Concentration       | Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.                                                                                                                                                                                                        |
| Conjugation         | Unconjugated                                                                                                                                                                                                                                              |
| Purification        | Affinity purified from tissue culture.                                                                                                                                                                                                                    |
| Dilution Range      | Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Immunofluorescence (0, 5-1, 0 µg/ml). Immunohistology (formalin-fixed: 1-2 µg/ml for 30 min at RT; staining of formalin-fixed tissues requires boiling tissue sections in 10mM citrate buffer, pH |
| Formulation         | PBS with 0.02% Sodium Azide.                                                                                                                                                                                                                              |
| Isotype             | IgG2ak                                                                                                                                                                                                                                                    |
| Storage Instruction | Store for up to 1 year at 2-8°C upon receipt.                                                                                                                                                                                                             |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|