

Anti-CD54-ICAM-1 antibody [EBS-CD-031] (STJ16101887)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	FC/FUNC/WB/IHC
Host / Source	Mouse
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Monoclonal
Concentration	Can be provided as 100 µg/mL, 500 µg/mL or 1mg/mL.
Conjugation	Unconjugated
Purification	Affinity purified from tissue culture.
Dilution Range	Flow cytometry (0, 5-1, 0 µg/million cells in 0, 1 ml). Functional studies (0, 02-2, 0 µg/ml without azide). Immunoblotting (0, 5-1, 0 µg/ml). Immunohistology (2-4 µg/ml for 30 min at RT; staining of formalin-fixed tissue requires boiling tissue sect
Formulation	PBS with 0.02% Sodium Azide.
Isotype	IgG2bk
Storage Instruction	Store for up to 1 year at 2-8°C upon receipt.
TARGET INFORMATION	
Gene ID	3383
Gene Symbol	ICAM1
UniProt ID	ICAM1_HUMAN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.