

Anti-CD44 antibody [F10-44-2] (STJ16101914)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IHC-F/FC/IF/IHC-P
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Purified
Formulation	PBS with <0.1% sodium azide
Isotype	IgG2a
Storage Instruction	Store aliquots at -20°C, and thawed aliquots at 2-8°C.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
UniProt ID	CD44_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---